Exhibit 25

**Boston
Scientific**

Guidezilla™

Guide Extension Catheter (5-in-6)

Directions for Use	2
Instrucciones de uso	4
Mode d'emploi	6
Gebrauchsanweisung	8
Istruzioni per l'uso	11
Gebruiksaanwijzing	13
Instruções de Utilização	16

GuidezillaTM

Guide Extension Catheter (5-in-6)

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your Boston Scientific representative.

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The Guidezilla guide extension catheter is a single lumen rapid exchange catheter compatible with 6F guide catheters and may be placed over an exchanged length 180 cm guidewire. The 145 cm device has a stainless steel proximal shaft with a 25 cm single lumen distal guide segment with a hydrophilic coating.

The Guidezilla device has two platinum-iridium marker bands, which enable visibility while using standard fluoroscopic methods. The distal marker band is located 2 mm from the distal tip. The proximal marker band is located 3 mm distal to the opening of the guide segment. The device has two positioning marks located at 90 cm (single mark) and 100 cm (double mark) from the distal tip, respectively.

The Guidezilla device is delivered through a guide catheter resulting in an inner diameter that is approximately 1 French size smaller than the guide catheter. The Guidezilla device has a proximal tab which indicates inner diameter and guide catheter compatibility.

Table 1. Compatibility Information

Model	Compatible Guidewire	Compatible Guide Catheter	Guidezilla Min. I.D.
6F (5-in-6)	≤ 0.014 in (0.36 mm)	≥ 6F / ≥ 0.070 in I.D. (1.78 mm)	.057 in (1.45 mm)

CONTENTS

Qty	Material
One (1)	Single-use Guidezilla guide extension catheter

INTENDED USE/INDICATIONS FOR USE

The Guidezilla guide extension catheter is intended to be used in conjunction with guide catheters to access discrete regions of the coronary and/or peripheral vasculature, and to facilitate placement of interventional devices.

CONTRAINDICATIONS

- Vessels less than 2.5 mm in diameter.
- Vessels in the neurovasculature and the venous system.

WARNINGS

- Use prior to the "use by" date as indicated on the label.
- Never advance the Guidezilla device into a vessel without a leading guidewire or without confirming location using fluoroscopic guidance. Vessel dissection or perforation may result.
- Due to the size and non-tapered tip of the Guidezilla device, extreme care must be taken to avoid vessel ischemia or vascular damage.
- Where there is limited clearance between interventional devices and the distal guide segment lumen, those devices must be advanced and withdrawn slowly with the hemostasis valve open to reduce the risk of embolism.
- This is a non-torqueable device. Torqueing the device may result in wire wrap or damage to the device or vessel.
- Never advance the Guidezilla device into a vessel with an effective diameter less than 2.5 mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the Guidezilla catheter, withdraw the Guidezilla catheter until the pressure returns to normal.

PRECAUTIONS

- Inspect the device prior to use for any bends or kinks. Any device damage may decrease the desired performance characteristics.
- The device lumen should be thoroughly flushed with heparinized saline prior to use.
- This device should be used only by physicians thoroughly trained in percutaneous, intravascular techniques and procedures.
- Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system. Use of systemic heparinization and heparinized sterile solution should be considered.
- Exercise care in handling of the device during a procedure to reduce the possibility of accidental breakage, bending or kinking.
- When the device is in the body, it should be manipulated only under fluoroscopy. Do not attempt to move the device without observing the resultant tip response.
- Never advance the Guidezilla device more than 15 cm beyond the tip of the guide catheter. Further distal advancement of the Guidezilla device could cause the entire guide segment to track outside of the guide catheter and impede withdrawal of the device.
- If strong resistance is encountered during manipulation of the devices, do not force passage. Determine the cause of the resistance before proceeding. If the cause cannot be removed, withdraw all devices simultaneously.

ADVERSE EVENTS

The adverse events include, but are not limited to:

- Vessel trauma (e.g., perforation, dissection)
- Vascular complication (e.g., puncture site complication)
- Thrombus
- Renal failure
- Slow-flow/occlusion
- Allergic reaction
- Death
- Emboli
- Hemorrhage/hematoma (e.g., vascular access complication)
- Infection
- Myocardial infarction
- Arterial spasm
- Intimal disruption

HOW SUPPLIED

Packaging is designed to maintain sterility according to expiration date on the label. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

Other items required but not provided:

- Guide catheter with an inner diameter large enough to accommodate the Guidezilla™ device
- Y-adaptor with hemostasis valve
- Guidewire with diameter ≤ 0.014 in (0.36 mm)
- Sterile syringe (for system flushing)
- Sterile heparinized saline (for system flushing)

Handling and Storage

Store in a cool, dry, dark place.

OPERATIONAL INSTRUCTIONS

Preparations for Use

1. Prior to use, carefully inspect the Guidezilla device packaging and components for damage.
2. Using sterile technique, transfer the protective tubing with the Guidezilla device into the sterile field.
3. Carefully remove the Guidezilla device from its protective tubing. Do not bend or kink the device during removal.
4. Immerse the distal guide segment in heparinized saline solution.

Delivery Procedure

Deliver the Guidezilla device according to the following steps:

1. Secure the previously inserted guidewire and backload the distal tip of the Guidezilla device onto the guidewire and advance until the device is just proximal to the hemostasis valve.
2. Open the hemostasis valve and advance the Guidezilla device through the hemostasis valve and into the guide catheter.

3. Under fluoroscopy, advance the Guidezilla device up to a maximum of 15 cm beyond the distal tip of the guide catheter and into the desired location within the vessel.

Warning: Due to the size and non-tapered tip of the Guidezilla device, extreme care must be taken to avoid vessel occlusion and damage to the wall of the vessels through which this catheter passes.

Warning: Never advance the Guidezilla device into a vessel with an effective diameter less than 2.5 mm. Vessel injury, ischemia, and/or occlusion may result. If pressure in a vessel dampens after inserting the Guidezilla catheter, withdraw the Guidezilla catheter until the pressure returns to normal.

4. Using fluoroscopy, confirm the desired position of the Guidezilla device in the vessel.
5. If performing an interventional procedure, backload the interventional device over the guidewire and advance the device through the guide catheter and Guidezilla device into the desired vascular space.

Note: Use caution when advancing the interventional device into the distal guide segment.

6. Tighten the Y-adaptor hemostasis valve securely on the proximal shaft of the Guidezilla device to prevent back-bleeding.
7. Perform the catheterization procedure. After completing the procedure, remove the Guidezilla device prior to removing the guide catheter from the vessel.

WARRANTY

Boston Scientific Corporation (BSC) warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond BSC's control directly affect the instrument and the results obtained from its use. BSC's obligation under this warranty is limited to the repair or replacement of this instrument and BSC shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. BSC neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **BSC assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

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© 2013 Boston Scientific Corporation or its affiliates. All rights reserved.
Boston Scientific International S.A.
Boston Scientific (Australia) Pty Ltd