
INSTRUCTIONS FOR USE
ROTAFLOW CENTRIFUGAL PUMP

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Rx ONLY CAUTION: US Federal law restricts this device to sale by or on the order of a physician.

CE 0124

TABLE OF CONTENTS

English		4
Deutsch		8
Français		12
Italiano		16
Español		20
Português		24
Polski		28
简体中		32
Performance Data		36

ROTAFLOW CENTRIFUGAL PUMP

1 DESCRIPTION

The ROTAFLOW Centrifugal Pump (RF-32) has a spinning rotor with flow channels which imparts rotary motion to the incoming blood, directing it through a spiral housing to the outflow port.

2 INDICATIONS

The RF-32 centrifugal pump is intended for use in an extracorporeal circulation during a cardiopulmonary bypass within the framework of surgical procedures, or in circulatory support. The RF-32 pumps blood within the flow rate range as described.

3 CONTRAINDICATIONS

Possible side effects include, but are not limited to: mechanical failure, hemolysis and thromboembolic phenomena. These are potential side effects with all extracorporeal blood systems. The centrifugal pump is contraindicated to use as a suction device.

4 WARNINGS AND PRECAUTIONS

- Only use the device in accordance with these Instructions for Use. Also follow the Instructions for Use provided for the console in use.
- Possible side effects include infections, hemolysis and embolisms. Side effects associated with extracorporeal circulation in general may also occur (e.g., post-perfusion syndrome and organ damage).
- Observe the use-by date on the device packaging.
- This disposable must not be reused or resterilized. Reused, reprocessed and resterilized devices could be contaminated, unsuitably packaged or unsuitably labeled. The material properties, function and sterility of the device may be compromised, resulting in failure, malfunction, incorrect storage or use of the device, non-compliance with the use-by date or infection of the patient or user. Damage to health of, or even death of, the patient, user or third parties may result.
- Do not use the device if it or the sterile packaging is damaged.
- Check system for completeness.
- Always observe strict asepsis when handling the device.
- Avoid contact of anaesthetic fluids (e.g. isoflurane, ethrane etc.) and alcoholic solutions with any surface of the centrifugal pump. Such contact could cause damage to the product. The operational safety is not guaranteed anymore and the manufacture's liability is out of the question.
- Insure that the pump and circuit has been debubbled and primed properly prior to beginning support.
- Attach tubing in such a manner as to prevent kinks or any restrictions that may alter flow.

- Make sure the complete system is watertight.
- Do not hit or strike the ROTAFLOW Centrifugal Pump with hard materials.
- To prevent a backflow (drainage) of patient blood, the pump must operate at a minimum speed or the tube must be clamped. The minimum speed is that at which the patient and line pressure are overcome.
- Please note the following when clamping the tube:
 - First, deactivate the flow regulation (⇒ ROTAFLOW System user's manual, section "Flow Regulation").
 - Reduce the pump speed to the minimum speed before clamping the tube, then turn the flow regulator to zero to prevent hemolysis.
 - When the pump is stopped, both the venous and arterial lines should be clamped. Inspect the arterial line for air bubbles that may arise as the result of retrograde flow. Remove any air bubbles detected.
- To restart the pump, again set it for the minimum speed so that both the patient and line pressure are overcome. Then remove the tube clamps. Never remove the clamps when the pump is operating at high speed otherwise the pump may magnetically decouple.
- While running the arterial line has to be monitored because the centrifugal pump, like all other centrifugal pumps, is able to pump air. The perfusion has to be stopped if air is detected. Gas embolism may cause severe health injuries.
- Appropriate heparinization is required and must be observed.
- Use an arterial filter.
- It is intended that systemic anticoagulation is utilized while this device is in use. Anticoagulation levels should be determined by the physician based on risks and benefits to the patient.
- If the outlet tubing is clamped the pump should not run for a long time to prevent the blood from being damaged.
- Frequent patient and device monitoring is recommended.
- Apply ultrasonic grease at both sides of the pump outlet to ensure ideal functioning of the flow measuring system and the bubble detector.
- If there are abnormal noises or blood is leaking the pump has to be changed.
- The adapter has to be cleaned and disinfected with a alcoholic moistened cloth, periodically.
- Dispose of the device according to the applicable regulations for contaminated devices and the packaging as per the national regulations.

If you have any questions or problems, or if you experience any malfunctions, please contact your local MAQUET representative or contact MAQUET Cardiopulmonary AG directly.

5 PACKAGING AND STORAGE

The device is delivered sterile and pyrogen-free. Sterility remains assured as long as the packaging is not opened or damaged and the use-by date has not expired. The device must be stored in a cool, dark and dry place until used.

Prescribed storage temperature: +10°C ... +30°C (+50°F ... +86°F)

6 SPECIFICATIONS

Priming Volume	approx. 32 ml
Inlet/Outlet I.D.	9.5 mm (3/8")
Max. rated pressure	750 mmHg (⇒ "Flow Pressure diagram", page 36)
Flow rates	0 ... 10 l/min

Materials

Housing, rotor, connector	Polycarbonate
Callote	Polyethylene
Ball	Sapphire
Magnet	Neodymium
Inlet and outlet caps 3/8"	HDPE

Adapter for the centrifugal pump RF-32

Model Number: MA-32 for Bio-Console from manufacture Medtronic.

7 DURATION OF APPLICATION

The ROTAFLOW Centrifugal Pump must be replaced after 4 days at the latest. Assess patient condition before replacing pump.

IMPORTANT: Check the ROTAFLOW Centrifugal Pump frequently for changes such as signs of coagulation, cracks and leaks. Perform a continuous visual inspection by illuminating the ROTAFLOW with a lamp. If in doubt, replace the centrifugal pump.

8 INSTALLATION

When using the ROTAFLOW Centrifugal Pump in the PLS Permanent Life Support System, it must not be driven by the Bio-Console (Medtronic) using the adapter MA-32. To run the RF-32 in the PLS System, always use only the ROTAFLOW Console RFC 20-970 with the corresponding RFD 20-973 Drive Unit.


The ROTAFLOW Centrifugal Pump is for use exclusively with the following systems.

- Centrifugal pump:
 - ROTAFLOW RF-32
- Centrifugal pump consoles:
 - ROTAFLOW Console RFC 20-970
 - Bio-Console (Medtronic) in combination with adapter MA-32 and the RF-32

As an alternative, the centrifugal pump console Bio-Console from Medtronic may, assuming it has 93/42 EEC approval, be used in conjunction with the appropriate magnetic adapter (MA-32).

Measuring systems when using Bio-Console (Medtronic):

- FLOWPROBE FP-32E
- Transducer TX 50 (Medtronic)

- 1 Mount the adapter as shown in figure 1.
- 2 After pressing the metal latch fix the pump onto the adapter.
- 3 Make sure the rim of the pump is underneath the metal rim of the adapter.

IMPORTANT: The user information of all products being used should be read carefully.

Make sure the outer edge of the ROTAFLOW is inserted below the guide nose of the adapter. The lid of the flow sample can still be closed if it is inserted incorrectly. However, the ROTAFLOW will then emit abnormal noises.

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
GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. ArjoHuntleigh focuses on patient mobility and wound management solutions. GETINGE provides solutions for infection control within healthcare and contamination prevention within life sciences. MAQUET specializes in solutions, therapies and products for surgical interventions and intensive care.

© MAQUET is a registered trademark of MAQUET GmbH & Co. KG - Order No. 70.000.42.43 - G.018 - 2010-06
Copyright by MAQUET Cardiopulmonary AG, Hirrlingen - Subject to modifications due to technical development