# Medtronic

## Bio-Pump™ BP-50, CBBP-50

Centrifugal Blood Pump; with or without Cortiva™ BioActive Surface

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Instructions for Use • Instrucciones de uso • Mode d'emploi

**! USA**

**Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

**CE**  
**0123**

The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.

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Achieve, Bio-Console, Bio-Pump, Cortiva, Medtronic

## Explanation of symbols • Explicación de los símbolos • Explication des symboles

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

## **1 Description**

The Bio-Pump centrifugal blood pump with or without Cortiva bioactive surface has been designed to utilize the constrained forced-vortex pumping principle. In action, blood is allowed to pass through a vortex created by the rotation of a series of smooth-surfaced rotating cones. Energy is transferred from the cones to the blood in the form of pressure and velocity as the blood is gently accelerated toward the outlet of the blood pump. The Bio-Pump blood pump is fabricated from sturdy, thromboresistant, biocompatible acrylics which are mounted on precision bearings for dependable operation. The non-occlusive hemodynamic design of the blood pump promotes laminar flow, allowing for improved blood handling capabilities and decreasing trauma which may be associated with extracorporeal circulatory support during cardiopulmonary bypass.

The blood pump couples to a safety-locked magnetic drive unit called the external drive motor (Model 540T), which is controlled by a Medtronic Bio-Console pump speed controller. The external drive motor and Bio-Console pump speed controller serve as the hardware components of the Medtronic extracorporeal blood delivery system.

The blood pumps have been sterilized using ethylene oxide.

Products coated with Cortiva bioactive surface include a “CB” prefix in the model number. Additional information regarding Cortiva bioactive surface is available in Chapter 7.

## **2 Indications for use**

The Medtronic Bio-Pump centrifugal blood pump is indicated for use only with the Medtronic Bio-Console pump speed controller to pump blood through the extracorporeal bypass circuit for extracorporeal circulatory support for periods appropriate to cardiopulmonary bypass procedures (up to 6 hours).

## **3 Contraindications**

The Medtronic centrifugal blood pumping system is contraindicated as a cardiotomy suction device.

Use the device only as indicated.

## **4 Adverse effects**

The following known adverse effects are associated with the use of the product: blood loss, coagulopathy, death, embolism, excessive blood component activation or thrombogenicity, exsanguination, hemolysis,

hemolytic anemia, infection, ischemia, neurological dysfunction, and organ dysfunction.

## 5 Warnings

Read all warnings, precautions, and instructions for use carefully before use. Failure to read and follow all instructions or failure to observe all stated warnings could cause serious injury or death to the patient.

- The blood pump has not been qualified through in vitro, in vivo, or clinical studies for long-term use (longer than 6 hours) as a bridge-to-transplant or for pending recovery of the natural heart.
- Possible side effects include, but are not limited to, infections, mechanical failure, hemolysis, and thromboembolic phenomena. These are potential side effects with all extracorporeal blood systems.
- Ensure that the blood pump and the circuit have been debubbled and primed properly before beginning bypass to minimize the risk of air reaching the patient. The use of an arterial filter is recommended.
- Massive air entry into the centrifugal blood pump deprimes the pump and stops blood flow. Stop the pump and remove air before resuming circulation.
- Do not use alcohol or alcohol-based solutions on any surface of the blood pump.
- To prevent retrograde blood flow when the blood pump outlet tubing is open, establish and maintain a minimum pump speed that overcomes line and patient resistance. Failure to do this could allow retrograde flow and exsanguinate the patient.
- Clamp the arterial lines when the blood pump is at rest to prevent retrograde flow.
- The blood pump is designed to be operated only with the Medtronic Bio-Console pump speed controller. There is no safety or performance data known to Medtronic which establishes compatibility of any other manufacturer's devices or components to the Medtronic system. Run the blood pump only on a properly maintained Bio-Console pump speed controller.
- Evaluate potential risk to the patient before changing the blood pump.
- Frequent patient and device monitoring is recommended; do not leave the centrifugal blood pump unattended while in use. Monitor the flow rate carefully for signs of occlusion throughout the circuit.
- Do not operate the centrifugal blood pump for more than 30 s in the absence of flow. The temperature within the pump could rise and increased cellular damage may result.

- The centrifugal blood pump must be handled using aseptic technique.
- Do not operate the blood pump with its inlet clamped. This will generate negative pressure in the blood pump and could cause air bubbles to form in the blood.

## **6 Precautions**

- This device should be used only by clinicians thoroughly trained in extracorporeal circulation procedures.
- Do not forcefully strike or drop the centrifugal blood pump. Shock may cause damage to the device, which may cause the device to malfunction.
- The fluid pathway of the centrifugal blood pump is sterile and nonpyrogenic. Inspect the device and package carefully before use. Do not use if the unit package or the product has been damaged or soiled, or if caps are not in place.
- This device was designed for single patient use only. Do not reuse, reprocess, or resterilize this product. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and create a risk of contamination of the device, which could result in patient injury, illness, or death.
- Attach tubing in a manner that prevents kinks or restrictions that may alter flow.
- Do not operate the centrifugal blood pump unprimed; damage to the internal components will occur.
- Do not use the blood pump if dropped. Shock may cause damage to the device, which may cause device malfunction.
- Do not use excessive force to install tubing on the pump because damage to the pump might occur.
- Take care to prevent damage to connectors when repositioning or attaching tubing.
- Always have a spare blood pump near and backup equipment available with appropriate protocol for change out.
- Do not place the centrifugal blood pump near items adversely affected by magnetic fields.
- Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.
- It is the user's responsibility to dispose of the devices in accordance with local regulations and hospital procedures.

## 7 Information for products with Cortiva bioactive surface

The primary blood-contacting surfaces of the product are coated with Cortiva bioactive surface. This coated surface enhances blood compatibility and provides a blood-contacting surface that is thromboresistant. Cortiva bioactive surface contains nonleaching heparin derived from porcine intestinal mucosa<sup>1</sup>.

**Caution:** Product coated with Cortiva bioactive surface is intended for single use only. Resterilization may adversely affect the Cortiva bioactive surface.

**Caution:** Follow a strict anticoagulation protocol and routinely monitor anticoagulation during all procedures.

## 8 Recommended priming procedures

1. Using accepted aseptic technique, attach the appropriate tubing to the inlet and the outlet of the blood pump.

**Note:** Applying a sterile saline solution to the inlet and outlet connectors may aid tubing attachment.

2. Where practical, flush the circuit and the blood pump with CO<sub>2</sub> (from filter or reservoir bag).
3. Using gravity, fill the blood pump with priming solution to a point beyond the flow probe and clamp. Eliminate air from the device by “walking” the air out of the outlet.

**Caution:** Do not hit or strike the blood pump with hands or instruments. Shock may cause damage to the device, which may cause device malfunction.

**Warning:** Ensure that the blood pump and the circuit have been debubbled and primed properly before initiating bypass to minimize the risk of air reaching the patient. The use of an arterial filter is recommended.

**Warning:** Massive air entry into the blood pump will cause it to deprime and blood flow to stop. Stop the blood pump and remove the air before resuming circulation.

4. With the outlet clamped, turn the Bio-Console pump speed controller power “ON.”
5. Zero (gain and balance) the flow probe according to the Bio-Console Operator and Reference Manual.
6. With the outlet still clamped, turn the RPMs of the Bio-Console pump speed controller to maximum RPM for 30 s. This delivers maximum

---

<sup>1</sup> Nonleaching is defined as heparin at a level  $\leq 0.1$  IU/mL, as measured through clinically relevant extraction conditions.

pressure within the blood pump, which aids in ensuring its integrity. Observe the blood pump for leakage or other anomalies.

**Warning:** Do not operate the blood pump for more than 30 s in the absence of flow. The temperature within the blood pump will rise and increased cellular damage may result.

7. Bring the RPMs to zero and recheck the integrity of the blood pump.

**Warning:** If leaks or other anomalies are found, remove the blood pump and replace with a new, sterile blood pump, repeating the above steps to prime.

8. If no anomalies are noted, continue priming the circuit.

9. Verify all connections and the integrity and flow of the circuit before use.

**Warning:** Do not operate the blood pump with its inlet clamped. This will generate negative pressure in the blood pump and could cause air bubbles to form in the blood.

## 9 Prebypass checklist

Information from this blood pump bypass checklist can be used to revise, upgrade, or expand existing checklists appropriately.

### 9.1 Equipment assembly

- Mount the motor drive correctly relative to the venous reservoir.
- Ensure that all of the electrical connections are secure.
- Test the control module power and display.
- Check the date and integrity of the sterile blood pump and disposable probe packages.
- Ensure that the flow transducer and the disposable probe are sized properly.
- Assemble the perfusion circuit in a sterile manner.
- Allow sufficient tubing length for the standby pumping unit.
- Connect the flow transducer and the disposable probe to the circuit in the correct location and flow direction.

### 9.2 Prime blood pump and circuit

- CO<sub>2</sub> flush the blood pump and the circuit, if indicated; turn off the CO<sub>2</sub>.
- Gravity prime and debubble the blood pump and the perfusion circuit.
- Check the blood pump for leaks, irregular motion, and noise.
- Check the circuit for visible air.
- Ensure that all of the tubing connections are secure.

- Clamp the blood pump outlet line completely.
- Clamp the venous return line completely.

### **9.3 Operating parameters**

- Calibrate the transducers according to the manufacturer's instructions.
- Set and verify the low and high flow alarms.

### **9.4 Emergency backup equipment**

Ensure that the following equipment is available:

- Backup power
- Handcrank (Model 150)
- Spare blood pump and disposable probe
- Control module or roller pump

### **9.5 Perfusion**

- Achieve minimum pump flow before unclamping lines.
- Monitor the control module for messages and alarms.
- Monitor the perfusion circuit for visible air.
- Maintain minimum pump flow before clamping lines.

### **9.6 Cleanup**

- Turn off the power.
- Properly discard disposable components.
- Cover the magnetic motor drive.
- Clean the console and the transducer.

### **9.7 Check equipment**

- Inspect and verify that the equipment is operational.
- Maintain the indicated preventative maintenance schedule.
- Ensure that the batteries are charged to full capacity.

## 10 Specifications

---

<b>Priming volume (approximate)</b>	48 mL	<b>Maximum flow rate</b>	1.5 L/min
<b>Inlet/outlet I.D.</b>	6.35 mm (1/4 in)	<b>Maximum blood pump speed</b>	4500 rpm
<b>Maximum operating pressure</b>	900 mmHg (120 kPa)	<b>Compatible systems</b>	Bio-Console pump speed controllers (with external drive motor Model 540T and emergency handcrank Model 150)

---

### 10.1 Storage and transport environment

Maximum transit temperature:

- BP-50: 65°C (150°F)
- CBBP-50: 40°C (104°F)

Store the product at room temperature.

## 11 Additional information

The following information and data are available upon request: test methods used to determine coverage, leaching, and bioactivity.

## 12 Limited warranty<sup>2</sup>

The following LIMITED WARRANTY applies to United States customers only:

- A. This LIMITED WARRANTY is available solely to the original purchaser of the Bio-Pump centrifugal blood pump with or without Cortiva bioactive surface, Models BP-50 and CBBP-50 (hereinafter, referred to as “product”).

- (1) In the event that the Product fails to function within normal tolerances due to a defect in materials or workmanship within a period of 1 year, commencing with the delivery of the Product to the purchaser, Medtronic will at its sole option: (a) repair or replace any

---

<sup>2</sup> This LIMITED WARRANTY is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only in the United States. Areas outside of the United States should contact their local Medtronic representative for exact terms of the LIMITED WARRANTY.

defective part or parts of the Product; (b) issue a credit to the purchaser equal to the purchase price (as defined in Subsection A(2) below) against the purchase price of the replacement Product; or (c) provide a functionally comparable replacement Product at no charge.

- (2) As used herein, purchase price shall mean the lesser of (i) the net invoiced price of the original Product, (ii) the price of currently functional comparable equipment, or (iii) the price of replacement equipment.

B. To qualify for the repair, replacement or credit set forth in Section A above, the following conditions must be met:

- (1) The Product must be returned to Medtronic within 60 days after discovery of the defect;
- (2) The Product must not have been used other than in strict accordance with the Instructions for Use, must not have been used on more than a single patient, and must not have been altered or repaired in any way which, in the judgment of Medtronic, affects the stability and reliability of the Product. The Product must not have been subjected to misuse, abuse or accident; and
- (3) The Use-by date for the Product must not have passed.

C. This LIMITED WARRANTY is limited to its express terms. In particular:

- (1) Except as expressly provided by this LIMITED WARRANTY, MEDTRONIC IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL, OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE, OR MALFUNCTION OF THE PRODUCT, WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT, OR OTHERWISE.
- (2) This LIMITED WARRANTY is made only to the purchaser of the Equipment. AS TO ALL OTHERS, MEDTRONIC MAKES NO OTHER WARRANTY, EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE WHETHER ARISING FROM STATUTE, COMMON LAW, CUSTOM, OR OTHERWISE. NO EXPRESS OR IMPLIED WARRANTY TO THE PURCHASER SHALL EXTEND BEYOND THE PERIOD SPECIFIED IN A(1) ABOVE. THIS LIMITED WARRANTY SHALL BE THE EXCLUSIVE REMEDY AVAILABLE TO ANY PERSON.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held to be illegal, unenforceable, or in conflict with applicable law by a court of

competent jurisdiction, the validity of the remaining portions of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

- E. No person has any authority to bind Medtronic to any representation, condition, or warranty except this LIMITED WARRANTY.

### **13 Limited warranty<sup>3</sup>**

**Important: This limited warranty does not apply in Australia.**

The following LIMITED WARRANTY applies to customers outside the United States:

- A. This LIMITED WARRANTY is available only to the original purchaser of the Bio-Pump centrifugal blood pump with or without Cortiva bioactive surface, Models BP-50 and CBBP-50 (hereinafter referred to as “Product”). In the event that the Product fails to function to the specifications set forth in the Instructions for Use, Medtronic will issue a credit, equal to the original Product purchase price (but not to exceed the value of the replacement Product) against the purchase of any Medtronic replacement Product used for that patient. **THE WARNINGS CONTAINED IN THE PRODUCT LABELING ARE CONSIDERED AN INTEGRAL PART OF THIS LIMITED WARRANTY. CONTACT YOUR LOCAL MEDTRONIC REPRESENTATIVE TO OBTAIN INFORMATION ON HOW TO PROCESS A CLAIM UNDER THIS LIMITED WARRANTY.**
- B. To qualify for the LIMITED WARRANTY, these conditions must be met:
  - (1) The Product must be used prior to its Use-by date.
  - (2) The Product must be returned to Medtronic within 60 days after use and shall be the property of Medtronic.
  - (3) The Product must not have been used for any other patient, must not have been misused, altered, abused or subject to accident, and the Product must not have been altered or repaired in any way which, in the judgment of Medtronic, affects its stability and reliability.
- C. This LIMITED WARRANTY is limited to its express terms. In particular:
  - (1) In no event shall any replacement credit be granted where there is evidence of improper handling, improper implantation or material alteration of the replaced Product.

---

<sup>3</sup> This LIMITED WARRANTY is provided by Medtronic, 710 Medtronic Parkway, Minneapolis, MN 55432. It applies only outside the United States.

- (2) Medtronic is not responsible for any incidental or consequential damages based on any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise.
- D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid.

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© 2018 Medtronic
Medtronic B.V.
Medtronic of Canada Ltd