



# MG 600 | Massage Gun with hot & cold function

Fast cooling as well as heating possible

Product number: 88987

Manual and other documents

## Description Reviews

### Product information "MG 600 | Massage Gun with hot & cold function"

#### Hot & cold function thanks to Peltier technology

The Hot&Cold function of the medisana MG 600 massage gun enables quick and effective cooling or warming of stressed muscles. It can be used both before training, where it prevents muscle strain through massage and the "hot function", and after training with the "cold function", which relaxes the muscles and reduces the risk of injury. The massage with the "hot function" also counteracts muscular hyperacidity and inflammation.

#### Self-massage on various parts of the body

Experience the relaxing deep tissue massage to activate muscles and promote muscle regeneration with the medisana MG 600 massage gun, which can be used to massage various parts of the body such as the shoulder, chest and back muscles as well as the gluteal, thigh and calf muscles. Muscle fatigue and sore muscles are a thing of the past. Three powerful vibration levels ensure deep muscle stimulation during use.

#### Perfect for at home, when travelling or during training




The Massage Gun MG 600 has no fewer than seven massage attachments for different muscle groups: spherical, U-shaped, curved, flat, spade-shaped, conical and the special Hot&Cold attachment. The practical storage case is perfect for use at home or on the go, for example during training, which also supports the long battery life of around 2-3 hours in normal operation. Thanks to its ergonomic shape and the accent colour "Stone-Blue", the MG 600 is also a real eye-catcher.

**Technical data**

- Size: 15 × 14.7 × 5 cm; Weight: 765 g
- Power supply: Lithium-ion battery (11.1 V; 2400 mAh); output power: 25 W
- Volume: < 60 dB (at a distance of 50 cm)
- Charging time: 4,5-5 hours
- Vibration levels: 1,800 rpm, 2,000 rpm, 2,400 rpm.

**Manufacturer Information**

**Kunden kauften auch:**



**Company**

medisana  
Health Magazine  
International  
Philosophy  
Job board

**Service**

Online Help  
Center  
Manuals  
Service contacts  
Returns  
Shipping  
Payment  
methods  
Contact

**Products**

Health Control  
Therapy &  
Healthy Home  
Flexible Heat  
Massage Devices  
Personal Care  
Sale  
Accessories

**Payment**



**Shipping**




**Social Media**




- Charging time: 4.5-5 hours
- Vibration levels: 1,800 rpm, 2,000 rpm, 2,400 rpm.

## Manufacturer Information


### Kunden kauften auch:





**FM 883 | Shiatsu foot massager**



**VC 150 | Cupping massager**



**AC 850 | Cellulite massager**

Company	Service	Products	Payment	Shipping	Social Media
<a href="#">medisana</a> <a href="#">Health Magazine</a> <a href="#">International</a> <a href="#">Philosophy</a> <a href="#">Job board</a>	<a href="#">Online Help Center</a> <a href="#">Manuals</a> <a href="#">Service contacts</a> <a href="#">Returns</a> <a href="#">Shipping</a> <a href="#">Payment methods</a> <a href="#">Contact</a>	<a href="#">Health Control</a> <a href="#">Therapy &amp; Healthy Home</a> <a href="#">Flexible Heat</a> <a href="#">Massage Devices</a> <a href="#">Personal Care</a> <a href="#">Sale</a> <a href="#">Accessories</a>			<a href="#">f</a> <a href="#">📺</a> <a href="#">@</a> <a href="#">in</a> <a href="#">🎵</a>



\* All prices incl. VAT, plus shipping costs | Free shipping from 40€ (to DE)